



GENOMICS

Pitfalls of primer and probe design and synthesis

Mary Span

Senior Product Manager

Eurogentec S.A.



EUROGENTEC
EGT GROUP



GENOMICS

Steps in qPCR assay

- Set up experiment
 - statistical relevant # samples/experimental group
 - controls
- **Design primers and probes**
- RNA extraction
 - quality of RNA
- Reverse Transcription reaction
 - one step or two step reaction
- qPCR reaction
 - singleplex or multiplex
- Data analysis



EUROGENTEC
EGT GROUP



GENOMICS

PCR efficiency

- High PCR efficiency
 - high accuracy
 - high reproducibility
- PCR efficiency influenced by
 - length of amplicon
 - GC content of amplicon
 - secondary structures in primers, probes, amplicons
 - concentration reaction components
 - PCR inhibitors/PCR enhancers
 - quality RNA/cDNA



EUROGENTEC
EGT GROUP



GENOMICS

PCR efficiency

- Easiest way to determine PCR efficiency: standard curve with R^2 close to 1,00 and intercept close to -3,32

Correlation Coefficient: 0.999 Slope: -3.317 Intercept: 8.156 $Y = -3.317 X + 8.156$
PCR Efficiency: 100.2 %

□ Unknowns
◇ Standards



Exponential amplification = $10^{(-1/\text{slope})}$

Efficiency = $10^{(-1/\text{slope})} - 1$



EUROGENTEC
EGT GROUP



GENOMICS

PCR efficiency

- 100% PCR efficiency
 - every PCR cycle amount of DNA is doubled
 - 2x dilution curve ΔC_t of 1 between every dilution
 - 10x dilution curve ΔC_t of 3,2 between every dilution
- Variation coefficient (R^2)
 - indication how well data points lie on one straight line
 - low R^2 indication for pipetting mistakes, inaccurate way of working, diluting out inhibitory factors



EUROGENTEC
EGT GROUP



GENOMICS

Why do you need a good primer and probe design?

- Well-designed primers and probes are a prerequisite for successful RT qPCR in terms of
 - high PCR efficiencies
 - specific PCR products
 - no co-amplification of genomic DNA
 - no amplification of pseudogenes
 - most sensitive results



EUROGENTEC
EGT GROUP



GENOMICS

Remarks primer and probe design

- Try to make primers fit the probe and not the probe fit the primers
- Do not expect that the primers you used before for your normal PCR will be the ideal pair to which a probe can be designed
- Using a design software is not a 100% guarantee to get a good primer/probe set, but is a good tool to make your life easier
- Especially with SYBR[®] green I assays; try several primer sets
- It is not always possible to design a primer/probe set for a specific sequence
- The first suggestion in the list of Primer Express[®] is the shorted amplicon, not the best primers and probe




EUROGENTEC
EGT GROUP



GENOMICS

Remarks primer and probe design

- Check T_m using several softwares
- If T_m differs by more than 3°C check T_m experimentally



Source: ABI User Bulletin 6 ABI PRISM® Sequence Detection System



GENOMICS

Design guidelines for SYBR[®] green I assays

- **Primers**
 - length
 - 18-30 bases
 - GC content
 - 30-80% (ideally 40-60%)
 - T_m
 - 63-67°C (ideally 64°C), so that Tannealing is 58-62°C (ideally 59°C)
 - ΔT_m forward primer and reverse primer < 4°C
 - avoid mismatches between primers and target, especially towards the 3' end of the primer
 - avoid runs of identical nucleotides, especially of 3 or more Gs or Cs at the 3' end
 - avoid 3' end T (allows mismatching)
 - avoid complementarity within the primers to avoid hairpins (check using a software)



EUROGENTEC
EGT GROUP



GENOMICS

Design guidelines for SYBR[®] green I assays

- **Amplicon**
 - length
 - 80-150 bp
 - shorter amplicons will give higher PCR efficiencies
 - longer amplicons will give a higher ΔR_n as more SYBR[®] green I is incorporated
 - GC content
 - 30-80% (ideally 40-60%)
 - avoid secondary structures in the amplicon (check with Mfold: www.bioinfo.rpi.edu/applications/mfold/)



EUROGENTEC
EGT GROUP



GENOMICS

Design guidelines for 5' exonuclease assays

- **Probes**
 - length
 - 18-30 bases
 - Optimal: 20
 - lengths over 30 bases are possible, but it is recommended to position the quencher not at the 3' end, but internally 18-25 bases from the 5' end
 - GC content
 - 30-80%
 - T_m
 - T_m of the probe must be 8-10°C (8°C for genotyping, 10°C for expression profiling) higher than the T_m of the primers
 - select the strand that gives the probe more Cs than Gs
 - place probe as close as possible to primers without overlapping them
 - avoid mismatches between probe and target
 - avoid runs of identical nucleotides, especially of 4 or more Gs
 - avoid 5' end G (quenches the fluorophore)
 - for multiplex assays: for genotyping
 - position the polymorphism in the center of the probe
 - adjust the probe length so that both probes have the same T_m



EUROGENTEC
EGT GROUP



GENOMICS

Design guidelines for 5' exonuclease assays

- **Primers**
 - length
 - 18-30 bases
 - GC content
 - 30-80% (ideally 40-60%)
 - T_m
 - 63-67°C (ideally 64°C), so that Tannealing is 58-62°C (ideally 59°C)
 - ΔT_m forward primer and reverse primer < 4°C
 - avoid mismatches between primers and target, especially towards the 3' end of the primer
 - avoid runs of identical nucleotides, especially of 3 or more Gs or Cs at the 3' end
 - avoid 3' end T (allows mismatching)
 - avoid complementarity within the primers to avoid hairpins (check using a software)



EUROGENTEC
EGT GROUP



GENOMICS

Design guidelines for 5' exonuclease assays

- **Amplicon**
 - length
 - 80-120 bp optimal (up to 1000 bp possible with adjusted reaction times)
 - shorter amplicons will give higher PCR efficiencies and more efficient 5' nuclease reactions
 - GC content
 - 30-80% (ideally 40-60%)
 - avoid secondary structures in the amplicon (check with Mfold: www.bioinfo.rpi.edu/applications/mfold/)




EUROGENTEC
EGT GROUP



GENOMICS

Example design 18S rRNA

- Comparison between two different primer-probe sets for 18S rRNA using same reaction components and experimental conditions





GENOMICS

Guidelines for fluorophores and quenchers

- Fluorophore has to fit real-time thermocycler

Thermocycler	Dye1	Dye2	Dye3	Dye4	Dye5	Dye6	Dye7
GeneAmp® 5700	FAM						
ABI Prism® 7000	FAM	VIC/YY/JOE	NED/TAMRA	ROX			
ABI Prism® 7700	FAM	VIC/YY/JOE/TET	NED/TAMRA	ROX			
ABI Prism® 7900	FAM	VIC/YY/JOE/TET	NED/TAMRA	ROX			
ABI Prism® 7300	FAM	VIC/YY/JOE	NED/TAMRA	ROX			
ABI Prism® 7500	FAM	VIC/YY/JOE	NED/TAMRA/Cy3	ROX/TR	Cy5		
i-cycler IQ®	FAM	VIC/HEX/TET/Cy3/YY	Cy3/TAMRA	ROX/TR	Cy5		
Mx3000P®	FAM	TET/YY	HEX/JOE/VIC/YY	TAMRA	Cy3	TR/ROX	Cy5/Alexa 350
Mx4000®	FAM	TET/YY	HEX/JOE/VIC/YY	TAMRA	Cy3	TR/ROX	Cy5
Rotorgene 2000	FAM	TET/JOE/VIC/YY	ROX/TAMRA/Cy3/TR	Cy5			
Rotorgene 3000	FAM	TET/JOE/VIC/YY	MAXROX/Cy3/TR	Cy5			
DNA Engine Opticon® 1	FAM						
DNA Engine Opticon® 2	FAM	TET/HEX/VIC/YY/TAMRA					
Chromo 4	FAM	TET/JOE/VIC/YY	ROX/TR	Cy5			
Smartcycler® 1	FAM	TET/JOE/VIC/YY	TAMRA/Cy3/Alexa	ROX/TR			
Smartcycler® 2	FAM	TET/Cy3/YY	ROX/TR	Cy5			
Lightcycler®	FAM	LC Red 640/ROX	LC Red 705/Cy5				
Lightcycler® 2.0	FAM	HEX/VIC/YY	LC Red 610	LC Red 640	LC Red 670	LC Red 705	
Quanta®	FAM	TET/HEX/VIC/YY/TAMRA					




EUROGENTEC
EGT GROUP



GENOMICS

Guidelines for fluorophores and quenchers

- Quencher has to fit fluorophore







GENOMICS

Guidelines for fluorophores and quenchers

- Singleplex or multiplex - TAMRA or DDQ1?




EUROGENTEC
EGT GROUP



GENOMICS

How to avoid or detect co-amplification of genomic DNA?

- Intron spanning primers




EUROGENTEC
EGT GROUP



GENOMICS

How to avoid or detect co-amplification of genomic DNA?

- Intron flanking primers



* Can be detected via melt curve





GENOMICS

How to avoid or detect co-amplification of genomic DNA?

- DNase I treatment of samples with RNase free DNase (Vandesompele, 2002) especially in cases of single exon genes
- Positions of exons and introns can be found in NCBI LocusLink databases (www.ncbi.nlm.nih.gov/LocusLink/)



EUROGENTEC
EGT GROUP



GENOMICS

How to avoid co-amplification of other genes?

- Check if primers are unique and specific by submitting your primers to a BLAST search (www.ncbi.nlm.nih.gov/BLAST/)



EUROGENTEC
EGT GROUP



GENOMICS

How to avoid primer-dimers or primer-probe dimers?

- Avoid complementarity between the primers, especially at 2 or more bases at the 3' ends of the primers (check using a software)
- Avoid complementarity of the probe with either of the primers (check using a software)



EUROGENTEC
EGT GROUP



GENOMICS

Recommended softwares

- Design of primers and probes
 - any primer design software/Oligo[®] 6.0/Primer 3.0 (http://frodo.wi.mit.edu/cgi-bin/primer3/primer3_www.cgi)
 - Primer Express[®] (TaqMan[®] probes)
 - BeaconDesigner[®] (Molecular Beacons)
 - Scorpio (Scorpions[®] primers)
- Verification of design
 - Mfold
 - BLAST
- **A software is just a tool to help you, not a guarantee for the perfect design!**



EUROGENTEC
EGT GROUP



GENOMICS

Additional information

- 📖 Your one-stop-shop real-time PCR supplier
- 📖 Troubleshooting guide for RT qPCR and qPCR
- 📖 Frequently asked questions for RT qPCR and qPCR
- 📖 This documentation is also available on www.eurogentec.com



EUROGENTEC
EGT GROUP



GENOMICS

Disclaimer

SYBR® is registered trademark of Molecular Probes Inc.

TaqMan® is a registered trademark of Roche Molecular Systems Inc.

Scorpions® is registered trademark of DxS Ltd.

ABI PRISM® is a registered trademark of The Perkin-Elmer Corp

iCycler iQ® is a registered trademark of Bio-Rad Laboratories Inc.

Mx3000p® and Mx4000® are registered trademarks of Stratagene

DNA Engine Opticon® is a registered trademark of MJ Research Inc.

SmartCycler® is a registered trademark of Cepheid Inc.

Primer Express® is a registered trademark of Applied Biosystems

Oligo® Primer Analysis software is a registered trademark of Molecular Biology Insights Inc.

BeaconDesigner® is a registered trademark of Premier Biosoft Inc.



EUROGENTEC
EGT GROUP