

# MUTATION DETECTION

Fragment Analyzer™ Automated CE System

Speed up reverse genetics lab work  
by accelerating sample prep and analyzing larger fragments.

Accurately detecting natural or induced-point mutations with slab gels is a long, drawn-out, time-consuming process. Now there's something better and faster. Mutation detection times are cut in half because of the Fragment Analyzer's speed and streamlined methods used with it.


## Mutation Detection Process Comparison – Fragment Analyzer™ Process vs Conventional Process

### FRAGMENT ANALYZER™ AUTOMATED CE SYSTEM PROCESS



- Step 1: PCR and heteroduplex formation
- Step 2: Heteroduplex transfer
- Step 3: Heteroduplex digestion
- Step 4: Diluent addition
- Step 5: Electrophoresis

### TRADITIONAL PROCESS



- Step 1: PCR and heteroduplex formation
- Step 2: Heteroduplex digestion
- Step 3: EDTA addition
- Step 4: Sephadex preparation and DNA capture
- Step 5: Acrylamide gel and apparatus preparation

- Step 6: Acrylamide gel prerun
- Step 7: Sample denaturation
- Step 8: Sample loading
- Step 9: Electrophoresis



## Compare post-PCR Process Steps for the Fragment Analyzer™ to the Slower Methods Used in Manual Acrylamide Gels.

---

### Mutation Detection Kit Process

1. Add 2µL PCR product to 2µL enzyme solution
2. Spin 10 sec
3. Incubate @45°C for 45 minutes
4. Add 20 µL diluent buffer and place on ice
5. Run CE

### Current Conventional Process

1. Add 20µL Cel I cocktail to PCR products
  2. Spin 1 min
  3. Incubate @45°C for 15 minutes
  4. Add 5µL of EDTA
  5. Purify on Sephadex or EtOH precipitate
  6. Run electrophoresis
- 

### Features and Benefits

#### No clean up step needed:

Eliminates several steps of the traditional process, reduces overall time and potential sample loss.

#### Potential to reduce gDNA input amount:

Smaller PCR setup and high sensitivity means less input gDNA is required — saves precious DNA.

#### Fast electrophoresis run times (40 minutes):

Get more separations done per instrument per day.

#### Minimal labor (no pouring gels or cleaning plates):

Automated process significantly reduces time handling fragile glass plates and toxic chemicals

#### Analyze up to 16 gene copies:

Maximize throughput with optimal organism pooling.

#### Examine fragments up to 10,000 base pairs:

Exceeds size limitations of traditional slab gel methods. Current slab gel methods can only analyze 1,500 base pairs.

#### Eliminate use of labeled primer sets:

Saves time and cost for expensive labels. No signal loss over time.

#### Ability to identify multiple cuts in one gene:

Sensitive intercalating dye allows easy detection of multiple fragment cut sites.

#### Analytical software for fragment sizing and concentration:

Easy to use data analysis software eliminates manual screening of gel pictures. Aids in displaying and sizing cut fragments

---

### Advanced Analytical Technologies, Inc.

2711 South Loop Drive, Suite 4150  
Ames, IA 50010 USA  
Phone: +1-515-296-6600  
Fax: +1-515-294-7141  
E-mail: sales-fs@aati-us.com  
[www.aati-us.com](http://www.aati-us.com)

### Advanced Analytical Technologies, GmbH.

Im Neuenheimer Feld 583  
D-69120 Heidelberg  
Phone: ++49 6221 868058-20  
Fax: ++49 6221 868058-99  
E-mail: sales-fs@aati-de.com  
[www.aati-de.com](http://www.aati-de.com)

