

# Quantitation of Plant miRNAs by RT-PCR

Karl J Guegler<sup>1</sup>, Julie Nguyen<sup>1</sup>, Cheng Lu<sup>2</sup>, Blake Meyers<sup>2</sup>, Pamela Green<sup>2</sup>, Kelly K McDonald<sup>1</sup>, Caifu X Chen<sup>1</sup>

<sup>1</sup> Advanced Research & Technology, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, CA 94404

<sup>2</sup> Department of Plant and Soil Sciences, Delaware Biotechnology Institute, Delaware Technology Park, 15 Innovation Way, Newark, DE 19711

## ABSTRACT

We have developed a novel, highly sensitive real-time PCR method for quantitation of miRNAs based on TaqMan<sup>®</sup> reagent chemistry. Assays have been designed and validated for over 300 known miRNAs from various species, including 50 from *Arabidopsis thaliana*. The sensitivity of this method allows detection of miRNAs in as little as 0.03 ng total RNA or 100 cells and the detection of Arabidopsis miRNAs in as little as 0.07 ng of total RNA. Using the synthetic Arabidopsis miRNA, miR159a, we were able to accurately determine C<sub>T</sub> values over seven orders of magnitude (R<sup>2</sup> = 0.995). This method allows us to discriminate between two miRNAs that differ by as little as a single nucleotide as well as between mature miRNAs and their precursors. The presence of genomic DNA or non-specific RNA does not affect the accuracy of the miRNA quantitation. Here we present data using this novel method to determine abundance and tissue-specific expression patterns of Arabidopsis miRNAs.

## INTRODUCTION

MicroRNAs are endogenous RNAs of approximately 22 nucleotides that play important regulatory roles in animals and plants by targeting mRNAs for cleavage or translational repression<sup>1</sup>. More than 700 miRNAs have been identified across species. Their expression levels vary greatly among species and tissues<sup>2</sup>. Low abundant miRNAs have been difficult to detect based on current technologies such as cloning, Northern hybridization<sup>3</sup>, and the modified Invader<sup>®</sup> assay<sup>4</sup>. Here, we present a new, real-time, quantitation method termed looped-primer RT-PCR for accurate and sensitive detection of miRNAs as well as other non-coding RNA (ncRNA) molecules. We have previously demonstrated the ability to accurately quantify miRNAs in human and mouse total RNA, cell lysates and whole cells. Here we demonstrate the ability to quantitate 50 miRNAs from five Arabidopsis tissue samples. For several miRNAs, we show comparative data from massively parallel signature sequencing (MPSS) and Northern analysis.

## MATERIALS AND METHODS

**miRNA targets:** A total of 50 miRNA assays were designed against Arabidopsis miRNA sequences downloaded from the Sanger database<sup>5</sup>.

**Tissue RNA samples:** Total RNA was isolated from seedling, root, callus, flower, and *rdr2* mutant flower tissues using a standard Trizol extraction method (Invitrogen).


**RT-PCR:** The assay includes two steps, RT and PCR (Figure 1). RT reactions (performed in duplicate) containing RNA samples, looped-primers, 1X buffer, reverse transcriptase, and RNase inhibitor were incubated for 30 minutes each, at 16°C and at 42°C. Real-time PCR reactions based on TaqMan<sup>®</sup> reagent chemistry were performed in triplicate on an Applied Biosystems 7900HT Fast Real-Time PCR System. No template controls (NTCs) were run for each assay under the same conditions.

**Data analysis:** The level of miRNA expression was measured using C<sub>T</sub> (threshold cycle). The C<sub>T</sub> is the fractional cycle number at which the fluorescence passes the fixed threshold. In the genomic cycle of PCR, the C<sub>T</sub> is inversely proportional to the sample starting copy number. The ΔC<sub>T</sub> was calculated by subtracting the C<sub>T</sub> of a control (internal or external) from the C<sub>T</sub> of the miRNA of interest. The ΔΔC<sub>T</sub> was calculated by subtracting the ΔC<sub>T</sub> of the calibrator (a sample used as the basis for comparative results) with the ΔC<sub>T</sub> of the test sample. Fold change was generated using the equation 2<sup>-ΔΔC<sub>T</sub></sup>.

**Northern:** RNA gel blots of small RNAs were carried out as described by Lu et al. Poster #767.

**MPSS:** MPSS is a RNA profiling method that has recently been adapted for the identification and quantification of small RNAs (Lu et al., P767).

Figure 1. Assay Scheme



## RESULTS


Table 1. Single Base Discrimination of miRNAs Assays

miRNA assay	Synthetic miRNA target						Relative detection (%) <sup>a</sup>
	let-7a	let-7b	let-7c	let-7d	let-7e	let-7f	
let-7a	100.0	0.3	7.2	0.1	0.0	0.0	100.0
let-7b	0.0	100.0	3.3	0.0	0.0	0.0	
let-7c	0.1	4.4	100.0	0.1	0.0	0.0	
let-7d	0.0	0.0	0.0	100.0	0.0	0.0	
let-7e	0.0	0.0	0.0	0.0	100.0	0.0	



(A) Relative detection (%) calculated based on C<sub>T</sub> difference between perfectly matched and mismatched assays. (B) Nucleotide sequences of let-7a RNA and closely related let-7 variants. Nucleotide differences are in red and denoted by an arrow. The number of mismatches between variants are shown to the right of the sequences.

Figure 2. Quantitation of Synthetic Arabidopsis miR159a



Amplification plots of fluorescence signal versus cycle number of synthetic Arabidopsis miR159a without mouse RNA (A), in the presence of 150 ng of non-specific mouse RNA (B), and combined data demonstrating the high degree of overlap (C). Estimated synthetic miRNA input (miR159a) in the RT step, based on OD, was 70, 700, 7,000, 70,000, 700,000, 7M, and 70M copies.

Figure 3. Sensitivity of miRNA Assays in Arabidopsis Seedling






Figure 4. Expression of miRNAs in Arabidopsis Flower



Histogram of data from 50 miRNAs detected in total RNA isolated from flower. Expression data were binned in two C<sub>T</sub> intervals, with the miRNAs in highest abundance grouped together in the 18.0-19.9 C<sub>T</sub> interval and miRNAs in lowest abundance grouped in the >32.0 C<sub>T</sub> interval. The graph shows a slightly skewed distribution with the majority of miRNAs in the 20.0-25.9 C<sub>T</sub> range.

Figure 5. Expression Ratios of miRNAs in Flower and Seedling




Histogram of data from 50 miRNAs detected in total RNA isolated from flower and seedling. The ratio of abundance of flower to seedling was calculated by dividing the average C<sub>T</sub> of flower by the average C<sub>T</sub> of seedling for each miRNA. The ratios were then binned to demonstrate the miRNAs that were most abundant in flower (<0.90 or 10% or greater C<sub>T</sub> in flower) or most abundant in seedling (>1.10 or 10% or greater C<sub>T</sub> in seedling). The largest number of miRNAs (bar shown in yellow) demonstrated relatively equal abundance in both flower and seedling.

Table 2. Data from MPSS Confirming RT-PCR Assays

miRNA	Flower (F)	Seedling (S)
miRNA167a,b,d	237,567	298,695
miRNA172a,b,c,d	7,492	226
mi399a,d,e,f	0	88


MPSS data from small RNA libraries prepared from flower (F) and seedling (S) are shown. miRNA167 and miRNA172 are examples of high abundance miRNAs and miRNA399 is an example of a low abundance miRNA as determined by RT-PCR of flower in Figure 4, and exhibit consistent MPSS data above. miRNA172 is more abundant in flower in both MPSS above and RT-PCR (Figure 5). miRNA399 is more abundant in seedlings in both experiments.

Figure 6. Northern Blot Validation of RT-PCR Assay



(Above) RNA gel blot hybridization of miRNA oligonucleotide probes to low mw total RNA from seedlings (S), flower (F), and flower from the *rdr2* mutant (*rd2*). The abundance of miRNA167 is similar in seedlings and flower, whereas for miRNA172 RNA abundance is much higher in flowers, similar to the RT-PCR data in Figures 4 and 5.

Figure 7. Heat Map Displaying Regulation of All miRNAs



(Right) Heat map displaying regulation of miRNAs in flower, seedling, and leaves compared to callus. The average C<sub>T</sub> value for miR159b was used to normalize data across callus, flower, seedling, and leaves, generating a ΔC<sub>T</sub> for each miRNA. ΔC<sub>T</sub> and fold changes were calculated as described in Materials and Methods using callus as the calibrator. The heat map shows miRNAs that are up- (red) and down- (green) regulated in flower, seedling and leaves as compared to callus. miR163 and miR171 showed the greatest up-regulation (>60-fold) and miR397a and miR397b demonstrated the greatest down-regulation (>12-fold) of flower, seedling, and leaves as compared to callus.

For more product information, please visit [miRNA.appliedbiosystems.com](http://miRNA.appliedbiosystems.com)

## REFERENCES

- Bartel, D. 2004. Cell 116, 281-297
- Kim, J. et al. 2004. PNAS 101:360-365
- Lim, L. et al. 2003. Gene & Development 17:991-1008
- Allawi, H. T. et al. 2004. RNA 10:1153-1161
- The microRNA Registry. Griffiths-Jones S. NAR, 2004, 32, Database Issue, D109-D111
- Lu et al. P767, PAG 2004

## ACKNOWLEDGEMENTS

Work at the University of Delaware was funded by SGER grant 0439186 from the NSF to P.J.G. and B.C.M.

## TRADEMARKS/LICENSES

For Research Use Only. Not for use in diagnostic procedures. The PCR process and 5' nuclease process are covered by patents owned by Roche Molecular Systems, Inc. and F. Hoffmann-La Roche Ltd. Applied Biosystems is a registered trademark and AB (Design) is a trademark of Applied Biosystems Corporation or its subsidiaries in the US and/or certain other countries. TaqMan is a registered trademark of Roche Molecular Systems, Inc. Invader is a registered trademark of Third Wave Technologies, Inc. All other trademarks are the sole property of their respective owners. © 2004 Applied Biosystems. All rights reserved.