



Quantitation using Real-Time PCR

Jennifer McMahon PhD

 **corbett**
LIFE SCIENCE
www.corbettresearch.com

 **corbett**
RESEARCH
www.corbettresearch.com

What is the aim of the Study?

To determine differences in gene expression levels between a control group and an experimental group(s).

Experiment considerations

Appropriate selection of Controls

Needs to be determined statistically

What are the groups being investigated?

Patients, treatments, time points

What genes are being targeted?

Based on micro-arrays, literature, physiological studies

Do I need a normalizer gene?

Can you control for the variation in starting template amount?

Can you rely on your U.V spectrophotometric results?

If not you need a normalizer

3

Statistics

Are you running a valid number of samples per group?

Are you running enough biological replicates per sample?

Choice of statistical test

- Student t-test (comparing two groups)
- 2 way ANOVA (number of groups over a number of variables)
- Repeated measures ANOVA (time course)
- Correlation analysis (relationship)

Outlier detection (Grubb's test)

4

Reverse Transcription

Convert total RNA to cDNA or just messenger RNA (mRNA)

Use of Random hexamer (RNA total) or OligodT primers (mRNA)

Avoid using oligodT primers if the target is at the 5' end of the RNA

Can't use oligodT if using 16S because it doesn't have a poly A tail

Choice of Enzymes available

RNAse H +/-

some work better on low copy template

one-step or two step reverse transcriptase real-time PCR?

Reaction Setup

General PCR rules apply e.g. don't set up reactions in electrophoresis areas

UNG - Uracil-DNA Glycosylase can prevent product contamination (not appropriate if product is to be used downstream)

NTCs - use no template controls to monitor contamination issues

RT-ve controls - use controls without reverse transcriptase to check for presence of DNA if using cDNA. Consider designing primers spanning exon/exon boundaries so that DNA can't be amplified

Normalizer or “Housekeeping” Genes

Normalize for variation in sample amount

Samples may vary in RNA extraction efficiency, RNA quality, cDNA synthesis, RNA (or DNA) concentration, pipetting of template

MUST be unregulated under experimental conditions

If amount of RNA \uparrow then amount of normaliser must \uparrow

i.e. normalizer must not change in response to treatment

Preferable to have similar abundance of normalizer & GOI

e.g. rRNA (18S) is more stable than mRNA and highly abundant so may not correlate well with low copy number genes

No such thing as a universal normalizer

Different experiments may affect normalizer expression

Selecting a normalizer

Check the literature, compare prospective normalizers, some people use 2 normalizers for all analysis

Vandestompe, J., DePreter, K., Pattyn, F., Poppe, B., Van Roy, N., De Paepe, A., Speleman, F. (2002) Accurate normalization of real-time quantitative RT-PCR data by geometric averaging of multiple internal control genes, *Genome Biology* **3(7)**

Choosing a Calibrator

Calibrator will be used for comparisons

Need a starting point or “normal”

GOI may be up regulated or down regulated as a result of “treatments” compared to the “normal”

Allows comparison of samples from various “treatments” over a number of runs – calibrator is like an internal control for run variation

Example of calibrator

Tissue culture experiment – untreated cells v treated cells

Pairwise normal tissue v tumour tissue – either choose one sample as a “normal” or pool a little from all “normals” and use as calibrator

Need a lot of calibrator

Calibrator will be used in every run to allow comparisons between runs

What do I use to calculate efficiency?

Important that the target used to calculate efficiency is similar to the samples

e.g. plasmid may amplify better than extracted material because it is cleaner (no inhibitors) and more pure (no other sequences present) for cross reactions

Serially dilute the target, plot and calculate the slope

The software does this for you, see next slide

9

Efficiency of PCR

Linear regression model

Serial dilution of the template


Plot the C_T vs log of concentration of template

Calculate the formula for the line of best fit

Efficiency is related to the slope (m)

$$E = (10^{-1/m}) - 1$$

$$E = 1 = 100\%$$



10

Quantitation

Relative to a Standard Curve

1. Absolute standard curve - results are in numbers
2. Relative standard curve - results are as ratios (comparative)


Relative to a reference sample (calibrator)

- Results are in ratios (up or down relative to calibrator)
- Normalizer used to correct for amount of template added
- The most powerful and widely used method
- Several formulae available

11

Absolute Quantitation

Standard curve can normalize C_T s to input amount



12

Comparative Quantitation

Available in Rotor-Gene Software only

Amplification Plot based Efficiency calculation

Based on the fluorescence history of each reaction


Software uses a second derivative of the raw amplification data

Software determines the “take off” point of a reaction – no need to draw a threshold

The slope of the line from the take off point until exponential amplification stops is used to calculate the amplification efficiency

Values are out of 2 – 2 is doubling. Anything above 1.6 generally OK

13



14

Notes

Amplification Value (Efficiency) of each reaction calculated

Average Amplification value used for analysis

Make sure you are looking at one gene at a time in case amplification efficiencies vary between genes


Variation in the Average Amplification value must be minimal

Switch off all NTCs and samples that have failed

Good assay quality control

15

Using Comparative Quantitation



Click on comp quant tab

Switch off NTCs and any samples that haven't worked as this will affect amplification value and std deviation

16

No.	C	Name	Take Off	Amplification	Rep. Takeoff	Rep. Amp.	Rep. Conc.	Rep. Calibrat
A2		Calibrator 6	17.1	1.86	17.1	1.78	1.00E+00	Calibrator
A3		Calibrator 6	17.1	1.75				
A4		Calibrator 6	17.0	1.72				
A5		Sample 2	17.9	1.97	17.8	1.82	6.72E-01	
A6		Sample 2	17.6	1.60				
A7		Sample 2	17.8	1.89				
A8		Sample 3	17.7	1.84	17.4	1.83	8.44E-01	
B1		Sample 3	17.1	1.71				
B2		Sample 3	17.3	1.94				
B3		Calibrator 5	18.6	1.79	18.5	1.71	4.44E-01	
B4		Calibrator 5	18.4	1.63				
B5		Calibrator 5	18.5	1.71				
B6		Sample 5	17.7	1.65	17.6	1.68	7.39E-01	
B7		Sample 5	17.3	1.68				
B8		Sample 5	17.8	1.71				
C1		Sample 6	18.0	1.96	17.8	1.71	6.60E-01	
C2		Sample 6	17.7	1.53				
C3		Sample 6	17.7	1.62				
C4		Sample 7	17.4	1.69	17.6	1.80	7.53E-01	

Calibrator Replicate

(A2) Calibrator 6

Results

Average Amplification 1.76±0.12

Want amplification values >1.6, SD as low as possible

Choose calibrator from drop down menu at side (defaults to first sample)

Export "Rep Conc" column to excel

Repeat analysis and calibrator selection for HK

Divide GOI by HK to correct for variation in starting amount

Relative Quantitation

Two standard curve method

2 Standard Curve Method

Generate a standard curve

for the GOI and for the HKG separately

use different channels or different pages to separate the curves

Calculate concentration

read values for the HKG and GOI of samples from their respective standard curves, divide one by the other

Standard Curves required every run

Rotor-Gene software can do analysis within a run but not between runs


Relative Quantitation

Comparative Ct ($\Delta\Delta$ Ct Method)


$\Delta\Delta C_t$ formula - $R = 2^{-\Delta\Delta C_T}$

ΔC_t = change in C_t ,

Compares the C_t difference for calibrator (GOI minus HK) and the sample (GOI minus HK)



C_T HK 9.47; C_T GOI 20.84; ΔC_T sample 11.37



C_T HK13.1: C_T GOI 20.66; ΔC_T 7.56

Livak, J. K., Schmittgen, T. D. (2001) Analysis of relative gene expression data using real-time quantitative PCR and the $2^{-\Delta\Delta C_t}$ method, *Methods* **25** p402-408


21

Using $\Delta\Delta C_T$ method


run standard curves for each gene at beginning of experiment

compare efficiencies of PCRs

efficiencies must be equal for this method




standard curve HKG



standard curve GOI

22

Assay validation



23

Rules and Assumptions

Formula assumes 100% efficiency

2 = doubling every cycle

Amplification efficiencies of GOI and HKG must be near identical

validation must be performed before using this method

assume thereafter that efficiency is constant so don't run standard curves every run

consider revalidation on semi-regular basis or with a new batch of reagents

Available in the Rotor-Gene software within a run but not between runs

24

Relative Quantitation REST method (Pfaffl)

25

REST Method (Pfaffl)

REST = **R**elative **E**xpression **S**oftware **T**ool

Formula

$$R = \frac{\text{Efficiency GOI}^{\Delta CP(\text{calibrator-sample})}}{\text{Efficiency HK}^{\Delta CP(\text{calibrator-sample})}}$$

Pfaffl efficiency out of 2, according to the formula $E = 10^{(-1/-m)}$

Rotor-Gene efficiency out of 1, according to the formula $E = (10^{(-1/-m)}) - 1$

To convert Rotor-Gene efficiency to Pfaffl efficiency add 1


26

Using REST Method (Pfaffl)


run standard curves for each gene at beginning of experiment

use software to calculate efficiency for HKG and GOI

compare efficiencies of PCRs – efficiencies can differ



standard curve housekeeper



standard curve GOI

27

Rules and Assumptions

Formula works on true efficiency

doubling not required or assumed

Amplification efficiencies of GOI and HKG can be different

efficiencies must be calculated prior to the analysis

assume thereafter that efficiency is constant so don't run standard curves every run

consider revalidation on semi-regular basis or with a new batch of reagents

Not available in the Rotor-Gene software –use spreadsheet or REST-RG


<http://www.gene-quantification.de/download.html#rest-2005>

28

Relative Quantitation Comparative Quantitation

29

Comparison of 4 methods of analysis



30

Reading

- Livak, J. K., Schmittgen, T. D. (2001) Analysis of relative gene expression data using real-time quantitative PCR and the $2^{-\Delta\Delta C_t}$ method, *Methods* **25** p402-408
- Pfaffl, M. W. (2001) A new mathematical model for relative quantification in real-time RT-PCR, *Nucleic Acids Research*, **29** p2002-2007
- Pfaffl, M. W. (2002) Relative expression software tool (REST®) for group-wise comparison and statistical analysis of relative expression results in real-time PCR, *Nucleic Acid Research*, **30** p 2-10
- Ramakers, C Ruijter, J. M., Lekanne Deprez, R. H., Moorman, A. F. M (2003) Assumption-free analysis of quantitative real-time polymerase chain reaction (PCR) data, *Neuroscience Letters* **339** p62-66
- Vandestompele, J., DePreter, K., Pattyn, F., Poppe, B., Van Roy, N., De Paepe, A., Speleman, F. (2002) Accurate normalization of real-time quantitative RT-PCR data by geometric averaging of multiple internal control genes, *Genome Biology* **3(7)** p

Analyzing a real run file...first steps

Test various primer sets on two samples


- Test two samples (in case one sample poorly extracted)
- Confirm that the primers do amplify the product (should run on a gel after amplification)

Check for efficiency and specificity

- Use comp quant to look at efficiency of each reaction
- Do melt analysis to compare true product and primer-dimer


Preliminary analysis using comp quant

- Can get an idea of whether the genes are up or down regulated between the two samples




Both samples have amplified well for this gene so can use the data
Need to remove NTCs from the analysis (lanes 16 & 52, amplification value of 0) to get tighter standard deviation for amplification value

33




Standard deviation for amplification value for gene 498 is good 0.03
Take off points are close (<0.3 cycles) = replicates are close.
Choose ywd2 as calibrator (drop down menu) or "1"; ywd3 has 1.68 fold more expression relative to ywd2 (as shown in the "Rep Conc" box).

34



Gene 550 has amplified well, standard deviation good, replicates close.
 If choose ywd2 as as calibrator then ratio is 1:1.68, exactly the same as for the last gene.
 If gene 498 was the HKG and gene 550 the GOI then there would be no difference in expression of 550 between the samples


35



Gene 766 has amplified well, standard deviation good, replicates close but a bigger spread than genes 498 and 550 (0.7 cycle spread for 766 compared with 0-0.4 cycles spread for 498 and 550)
 If choose ywd2 as a calibrator then ratio is 1:2.44, different to the ratio with genes 498 and 550 (1:1.68) but not a significant difference due to cycle spread.

36


Poor amplification



Both samples have amplified very poorly for this gene so can't use the data

37

Unequal amplification between samples



Ywd2 has amplified far better than Ywd3 for this gene so can't use the data

38