

# New concepts for accelerated real-time PCR analysis



Thorsten Traeger,<sup>1</sup> Kirsten Hildebrand,<sup>1</sup> Dirk Loeffert,<sup>1</sup> Ralf Peist,<sup>1</sup> Annette Tietze,<sup>1</sup> Andrea Stutte-Rabe,<sup>1</sup> Ulla Deutsch,<sup>1</sup> Dile Holton,<sup>2</sup> and Andreas Missel<sup>1</sup>

<sup>1</sup> QIAGEN GmbH, Hilden, Germany; <sup>2</sup> QIAGEN Sciences, Germantown, MD, USA


## Introduction

For researchers needing to increase their throughput or share a cycler with other users, there is a strong demand for faster, real-time PCR.

**Fast, real-time PCR can be achieved by:**


- Reduced DNA polymerase activation time
- Shortened amplification cycles
- Combined annealing and extension steps
- Use of a dedicated fast-cycling instrument
- Shortened RT step in one-step RT-PCR

Until now, fast, real-time PCR using standard reaction chemistry has been hampered by reduced sensitivity and increased variability of quantification data (1). We demonstrate how the combination of a newly developed fast-cycling PCR buffer with a rapid-activating hot-start DNA polymerase allow significant reduction of PCR cycling times without sacrificing specificity and sensitivity.




## New chemistries for fast-cycling PCR

We have developed a fast-cycling PCR buffer that significantly reduces denaturation, annealing, and extension times. A novel additive, Q-Bond, dramatically increases the binding affinity of DNA polymerase to single-stranded DNA. This turns the 3-step process of template denaturation, primer annealing, and DNA polymerase binding in standard-cycling PCR (A) into a faster 2-step process (B).



High annealing specificity is maintained by a balanced combination of KCl and NH<sub>4</sub>Cl in the buffer. The binding of primers to imperfectly matching sequences on the template is suppressed.



A novel enzyme, HotStarTaq® Plus DNA Polymerase, is rapidly activated in 3 or 5 minutes by a 95°C incubation at the start of PCR.

## Ultrafast cycling for end-point PCR analysis


Although fast results in PCR can be achieved on cyclers with rapid ramping rates, even faster results are possible by reducing cycling times. The QIAGEN® Fast Cycling PCR Kit, which integrates the fast-cycling PCR buffer with HotStarTaq Plus DNA Polymerase, provides significant time savings of up to 78% in end-point PCR. Fast results can be accomplished on all cyclers, including cyclers not capable of fast ramping rates.

The PCR buffer minimizes amplification of nonspecific products, primer-dimer formation, and background smear in every PCR cycle. Q-Solution, an additive that enables efficient amplification of "difficult" (e.g., GC-rich) templates, is also provided with the kit.

**PCR cycling times calculated for different fragment lengths\***


Fragment length	QIAGEN Fast Cycling procedure (min)	Standard cycling procedure (min)	Time saving
200 bp	15	68	78%
500 bp	20	68	71%
1000 bp	29	85	66%
3000 bp	63	155	59%

\*Total time required for a PCR run of 35 cycles. The specified PCR cycling times do not include ramping times, which are cycler-dependent.




## Fast SYBR Green quantification with high specificity

In biological systems, minute changes in transcript abundance often lead to strong biological effects. Therefore, a method for reliable and reproducible discrimination between similar copy numbers is critical. With QuantFast™ SYBR Green Kits, even small differences in the amount of low-copy targets can be clearly distinguished.




## Speed and sensitivity in probe-based detection

Reducing the duration of each PCR cycle leads to faster PCR, but can impair PCR performance. The use of our fast-cycling technology by QuantFast Probe Kits enables faster results in probe-based real-time PCR without compromising sensitivity. Fast and sensitive quantification is possible on all available cyclers and with different types of sequence-specific probe, such as TaqMan® and FRET probes.



## Summary

- A patent-pending, fast-cycling PCR buffer containing Q-Bond significantly reduces denaturation, annealing, and extension times.
- HotStarTaq Plus DNA Polymerase possesses no enzyme activity prior to PCR, and is rapidly heat-activated in 3 or 5 minutes.
- The fast-cycling PCR conditions provide significant time savings of up to 78% in end-point PCR and up to 60% in real-time PCR without compromising specificity and sensitivity.
- New QIAGEN chemistries enable fast-cycling on all cyclers, including those not capable of achieving rapid ramping rates.



**References**  
1. Hilscher, C., Vahron, W., and Dittmer, D.P. (2005) Faster quantitative real-time PCR protocols may lose sensitivity and show increased variability. Nucleic Acids Res. 33:e182.

Trademarks: QIAGEN®, HotStarTaq®, QuantFast®, QIAGEN Group, ABI PRISM®, Applied Biosystems® (Applied Biosystems or its subsidiaries), LightCycler®, TaqMan® (Roche Group), Mastercycler® (Eppendorf AG), SYBR® (Molecular Probes, Inc.). Purchase of the product QuantFast SYBR Green PCR Kit and QuantFast SYBR Green RT-PCR Kit is accompanied by a limited, nontransferable license from us to use the product for a defined period of time as described in U.S. Patent Nos. 5,994,056 and 6,171,785 and corresponding patent claims outside the United States for the purchaser's own internal research. No resale, reproduction, or other patent rights or any other patent rights, and no right to use the product for any other purpose are granted separately, by implication or by assignment. Further information on purchasing licenses may be obtained by contacting the Director of Licensing, Applied Biosystems, 850 Lincoln Center Drive, Foster City, California 94044, USA. Purchase of this product QuantFast SYBR Green RT-PCR Kit is accompanied by a limited, nontransferable license under RT and Reverse Transcriptase/PCR patents owned by Roche Molecular Systems, Inc. and T. Hoffmann in Roche AG to use the product for any other purpose is hereby granted separately, by implication or by assignment. No right to use the product for any other purpose is hereby granted separately, by implication or by assignment. © QIAGEN, all rights reserved.