



## Instruction Manual

# **SuperScript™ III CellsDirect cDNA Synthesis System**

Catalog Nos. 18080-200 and 18080-300

**Version A**  
14 May 2004  
25-0731



Table of Contents

Kit Contents and Storage..... iv

Introduction.....v

Lysing Cells .....1

First-Strand cDNA Synthesis.....5

PCR and Quantitative PCR .....6

Troubleshooting..... 10

References ..... 11

Related Products ..... 12

Purchaser Notification ..... 13

Technical Service ..... 15

# Kit Contents and Storage

## Shipping and Storage

Kit components are shipped on dry ice and should be stored at -20°C.

## Kit Components

Catalog no. 18080-200 provides reagents for 25 reactions.  
Catalog no. 18080-300 provides reagents for 100 reactions.

<u>Component</u>	<u>25 Rxns</u>	<u>100 Rxns</u>
Resuspension Buffer	250 µl	1 ml
RNaseOUT™ Recombinant		
Ribonuclease Inhibitor (40 units/µl)	50 µl	200 µl
DNase I (1 U/µl)	125 µl	500 µl
10X DNase I Buffer	40 µl	160 µl
25 mM EDTA	30 µl	120 µl
Oligo(dT) <sub>20</sub> (50 µM)	50 µl	120 µl
10 mM dNTP Mix	25 µl	100 µl
SuperScript™ III RT (200 units/µl)	25 µl	100 µl
5X RT Buffer*	150 µl	600 µl
0.1 M DTT	50 µl	100 µl
<i>E. Coli</i> RNase H (2 U/µl)	30 µl	100 µl
HeLa Total RNA (10 ng/µl)	10 µl	10 µl
Forward Control Primer (10 µM)	10 µl	10 µl
Reverse Control Primer (10 µM)	10 µl	10 µl

\* 5X RT Buffer composition: 250 mM Tris-HCl (pH 8.3, room temp.), 375 mM KCl; 15 mM MgCl<sub>2</sub>

## Materials Supplied by the User

The following additional items are required for use with this kit:

- Coulter Counter or hemacytometer
- Microcentrifuge
- Thermal cycler
- Trypsin (for adherent cell cultures only)
- 1X cold phosphate-buffered saline PBS, without calcium or magnesium
- 0.2-ml thin-walled PCR tubes or 96-well PCR plates
- Ice
- Pipettes
- Disposable gloves



### Note

DNA polymerase is not included in this kit. Recommended DNA polymerases and optional amplification protocols are provided starting on page 6.

# Introduction

---

## System Overview


The SuperScript™ III CellsDirect cDNA Synthesis System is an optimized kit for synthesizing first-strand cDNA directly from mammalian cell lysate without first isolating the RNA. Lysis and reverse transcription are performed in the same tube, and the resulting first-strand cDNA is ready to use in cloning, PCR, or real-time quantitative PCR (qPCR) applications.

In traditional RT-PCR, RNA is first isolated from cells in a time-consuming procedure that can lead to a loss of material. Using the SuperScript™ III CellsDirect cDNA Synthesis System, the cells are lysed and the cDNA is generated from the lysate in a single tube with minimal handling and no sample loss. DNase I is added to eliminate genomic DNA prior to first-strand synthesis.

This kit has been optimized for small cell samples, ranging from 10,000 cells down to a single cell (as measured by serial dilution). The use of SuperScript™ III Reverse Transcriptase ensures high specificity and high yields of cDNA from small amounts of starting material—as little as 10 pg total RNA.

After synthesis, the first-strand cDNA can be amplified with specific primers by PCR without intermediate organic extractions or ethanol precipitations.

The diagram below outlines the procedure:



---

*Continued on next page*

## Introduction, Continued

---

### Advantages of the Kit

This kit has the following advantages:

- Compatible with a wide range of mammalian cell types grown under different treatment conditions
  - Single-tube format minimizes reagent loss, sample loss, and handling time
  - Total lysate volume is used in first-strand cDNA synthesis reaction, providing greater yields with a limited number of cells and allowing for detection of rare transcripts
  - SuperScript™ III Reverse Transcriptase, with reduced RNase H activity and higher thermal stability, produces high yields of cDNA in the first-strand synthesis reaction, for greater sensitivity and enhanced detection of rare transcripts
  - Generates high-quality cDNA for use in a variety of applications, including cloning, PCR, and real-time qPCR
  - Simple protocol takes less than 2 hours
- 

### SuperScript™ III RT

SuperScript™ III Reverse Transcriptase is an engineered version of M-MLV RT with reduced RNase H activity and increased thermal stability. The enzyme can be used to synthesize first-strand cDNA at temperatures up to 55°C, providing increased specificity, higher yields of cDNA, and more full-length product than other reverse transcriptases.

Because SuperScript™ III RT is not inhibited significantly by ribosomal and transfer RNA, it can effectively synthesize first-strand cDNA directly from total RNA. The concentration of SuperScript™ III RT in this system has been optimized to synthesize first-strand cDNA from total RNA in cell lysate.

---

### Control RNA and Primers

The control RNA provided with this system consists of HeLa Total RNA (10 ng/μl). The Forward Control Primer and Reverse Control Primer provided with this kit are designed from the human GAPDH gene and produce a 1.18-kb PCR product.

---

### Quality Control

This kit was verified in an end-point RT-PCR reaction using a serial dilution of HeLa cells and GAPDH primers, yielding a 1.18 kb PCR product. No template controls and negative reverse transcriptase controls were also prepared. Results were confirmed by gel electrophoresis.

---

# Methods

## Lysing Cells

### Introduction

In this step, you lyse your cells in Resuspension Buffer and perform a DNase I digestion to remove genomic DNA from the sample.

### Cell Types and Density

This kit has been optimized for small cell samples, ranging from 1 to 10,000 cells. The performance of this kit was verified using several different mammalian cell lines, including HeLa, COS-7, 293, Jurkat, CV1, and K562. Cells may be grown under variety of conditions and treatments. Any type of culture vessel can be used.



### Important

- We recommend using a maximum of 10,000 cells per reaction. Higher numbers of cells may inhibit reverse transcription and result in reduced yields and/or truncated cDNA product.
- Make sure that all solutions and equipment that come in contact with the cells are sterile. Always use proper sterile technique and work in a laminar flow hood when handling cells.

### Required Materials

The following materials are provided by the user:

- Mammalian cell cultures in growth media
- Coulter Counter or hemacytometer
- Centrifuge (for pelleting cells)
- Incubator, water bath, or thermal cycler preheated to 75°C
- Trypsin (for adherent cell cultures only)
- 1X cold phosphate-buffered saline (PBS), without calcium or magnesium
- 0.2-ml thin-walled PCR tubes or 96-well PCR plates
- Ice
- Pipettes

The following materials are provided in the kit:

- Resuspension Buffer
- RNaseOUT™ (40 U/μl)
- DNase I, Amplification Grade (1 U/μl)
- 10X DNase I Buffer
- EDTA, 25 mM
- Optional: Control HeLa Total RNA

### Control Reaction

For the control reaction, use 1 μl of the HeLa Total RNA provided in the kit instead of cell lysate.

*Continued on next page*

## Lysing Cells, continued

---



### Note

All steps should be performed on ice, and reagents should be chilled and/or thawed immediately prior to use.

The incubator should be **preheated** to 75°C.

---

### Lysing Adherent Cells or Cells in Suspension

Use the following lysis procedure for adherent cell cultures in vessels larger than 24-well plates. For cells in suspension, skip Steps 1–4 and proceed to Step 5 below.

1. Add enough trypsin to cover the adherent cells in your tissue culture dish, plate, or flask (e.g., for a 10-cm dish, use ~1 ml; for a T75 flask, use ~3 ml).
2. Incubate for 5 minutes at room temperature or in a 37°C incubator.
3. Check for cell detachment under a microscope. If cells have not detached, gently tap the dish or flask to dislodge the cells, or let the cells incubate longer, checking them every minute under a microscope.
4. When all the cells have detached, add serum-containing media to a final volume of 10 ml (for 6- and 12-well plates, add a 1X–2X volume of media). Note that the media must contain serum to inactivate the trypsin.
5. Pipet the cells gently up and down to mix, and then transfer the cell suspension to a centrifuge tube.
6. Spin the cells at  $200 \times g$  for 5 minutes to pellet.
7. Aspirate the media and wash the cell pellet with 5–10 ml of 1X cold PBS.
8. Spin the cells at  $200 \times g$  for 5 minutes to pellet.
9. Aspirate the PBS and resuspend the pellet in 500  $\mu$ l to 1 ml of 1X cold PBS. Mix the cell solution gently.
10. Collect a small aliquot to verify that the cells at the desired concentration. Determine cell density electronically using a Coulter Counter or manually using a hemacytometer chamber.
11. Adjust the cell density using cold PBS so that it falls within the range of 1–10,000 cells/ $\mu$ l. Count the cells again to verify cell concentration.
12. To a 0.2-ml thin-walled PCR tube or plate well **on ice**, add 1  $\mu$ l of RNaseOUT™ (40 U/ $\mu$ l) and 10  $\mu$ l of Resuspension Buffer.
13. Transfer 1–2  $\mu$ l of cells (<10,000 cells) to the PCR tube/well.  
**Control:** For the control reaction, add 1  $\mu$ l of the Control HeLa Total RNA to the PCR tube or plate well instead of cell lysate.
14. Transfer the tube/plate to an incubator, water bath, or thermal cycler preheated to 75°C and incubate for 10 minutes. **Control:** For the control reaction, incubate for 3 minutes.
15. After incubation, spin briefly to collect the condensation and proceed to **DNase I Digestion**, page 4.

---

*Continued on next page*



## Lysing Cells, continued

---

### Lysing Cells in Tissue Culture Wells

**Note:** Seed cells in tissue culture wells so that 10  $\mu$ l of resuspended cells will yield the desired concentration.

For adherent cells grown in tissue culture wells (i.e., in 24-well, 48-well, or 96-well plates), perform the following lysis procedure.

1. Aspirate the media in each well and wash each well with 1X cold PBS. Aspirate the PBS.
2. Add Resuspension Buffer to each well. For 96-well plates, add at least 10  $\mu$ l of buffer to each well. For 24-well plates, add at least 100  $\mu$ l of buffer to each well. The buffer should cover the cells in the well.
3. Incubate the plates on ice for up to 10 minutes. During that period, tap the plate periodically and check the cells under a microscope every 2–3 minutes to see whether they have detached or burst.
4. After 10 minutes, gently pipet the cells up and down to dislodge the remaining attached cells. Count the cells or estimate their density based on the seeding density (10  $\mu$ l should contain <10,000 cells).
5. Transfer 10  $\mu$ l of the cell suspension to a 0.2-ml thin-walled PCR tube or plate well.

**Control:** For the control reaction, add 10  $\mu$ l of Resuspension Buffer to a PCR tube or plate well, and then add 1  $\mu$ l of Control HeLa Total RNA.

6. Add 1  $\mu$ l of RNaseOUT™ (40 U/ $\mu$ l) to the PCR tube/well.
7. Transfer the tube/plate to an incubator or thermal cycler preheated to 75°C and incubate for 10 minutes. **Control:** For the control reaction, incubate for 3 minutes.
8. After incubation, spin briefly to collect the condensation, and proceed to **DNase I Digestion**, page 4.

---

*Continued on next page*

# Lysing Cells, continued

---

## DNase I Digestion

In this step, you treat the cell lysate with DNase I to degrade any contaminating DNA.

1. Place each tube/plate from Step 15, page 2, or Step 8, page 3, on ice, and add the following:

<u>Component</u>	<u>Amount</u>
DNase I, Amplification Grade (1 U/μl)	5 μl
10X DNase I Buffer	1.6 μl

2. Mix by gently pipetting up and down, and spin briefly to collect the contents.
  3. Incubate for 5 minutes at room temperature. **Note:** A longer incubation time (up to 10 minutes) may be used for larger samples (>1,000 cells). However, incubation times exceeding 10 minutes can greatly reduce cDNA yield.
  4. Spin briefly, and add 1.2 μl of 25 mM EDTA to each tube/well on ice. Mix by gently pipetting up and down, and spin briefly to collect the contents.
  5. Incubate at 70°C for 5 minutes.
  6. Spin briefly and proceed to **First-Strand cDNA Synthesis**, page 5.
-

# First-Strand cDNA Synthesis

---

## Required Materials

The following materials are provided by the user:

- Thermal cycler preheated to 70°C
- Ice
- Pipettes

The following materials are provided in the kit:

- Oligo(dT)<sub>20</sub> (50 µM)
  - 10 mM dNTP Mix
  - 5X RT Buffer
  - RNaseOUT™ (40 U/µl)
  - SuperScript™ III RT (200 U/µl)
  - 0.1 M DTT
  - RNase H (2 U/µl)
- 

## First-Strand cDNA Synthesis

1. Place each tube from **DNase I Digestion**, Step 6, page 4, on ice, and add the following:

<u>Component</u>	<u>Amount</u>
Oligo(dT) <sub>20</sub> (50 mM)	2 µl
10 mM dNTP Mix	1 µl

2. Mix by gently pipetting up and down, and spin the tube briefly to collect the contents.
3. Incubate the tube at 70°C for 5 minutes. Spin the tube briefly to collect the condensation.
4. Place the tube on ice for 2 minutes, and then add the following:

<u>Component</u>	<u>Amount</u>
5X RT Buffer	6 µl
RNaseOUT™ (40 U/µl)	1 µl
SuperScript™ III RT (200 U/µl)*	1 µl
0.1 mM DTT	1 µl

\*For negative RT controls, use 1 µl of sterile, distilled water instead of SuperScript™ III RT

5. Mix by gently pipetting up and down, and spin the tube briefly to collect the contents.
  6. Transfer the tube to a thermal cycler preheated to 50°C. Incubate for 50 minutes.
  7. Inactivate the reaction at 85°C for 5 minutes.
  8. Add 1 µl of RNase H (2 U/µl) to each tube and incubate at 37°C for 20 minutes. **Note:** This step is optional if you are amplifying short targets (<1.0 kb) in end-point PCR; it is **required** for qPCR.
  9. Chill the reaction on ice.
  10. Store the single-stranded cDNA at -20°C, or proceed directly to PCR amplification.
-

# PCR and Quantitative PCR

**Introduction**      The first-strand cDNA generated using this kit can be used directly in PCR or real-time quantitative PCR (qPCR) without additional purification. This section provides example protocols for PCR and high-fidelity PCR.

<b>PCR and Real-Time qPCR Kits</b>	The following PCR enzymes and qPCR kits are available separately:		
	<u>Product</u>	<u>Size</u>	<u>Cat. No.</u>
	Platinum® <i>Taq</i> DNA Polymerase	100 rxns	10966-018
		250 rxns	10966-026
		500 rxns	10966-034
		5000 rxns	10966-083
	Platinum® <i>Taq</i> DNA Polymerase High Fidelity	100 rxns	11304-011
		500 rxns	11304-029
		5000 rxns	11304-102
	Platinum® Quantitative PCR SuperMix-UDG	100 rxns	11730-017
		250 rxns	11730-025
	Platinum® SYBR® Green qPCR SuperMix-UDG	100 rxns	11733-038
		250 rxns	11733-046

**PCR Enzymes**      For amplifying the first-strand cDNA generated using this kit, we recommend Platinum® *Taq* DNA Polymerase for targets < 1.0 kb and Platinum® *Taq* DNA Polymerase High Fidelity for targets > 1.0 kb.

**Platinum® *Taq* DNA Polymerase** is recombinant *Taq* DNA polymerase complexed with proprietary Platinum® antibodies that block polymerase activity at ambient temperatures. Activity is restored after the denaturation step in PCR cycling at 94°C, providing an automatic “hot start” for *Taq* DNA polymerase in PCR. Hot starts in PCR provide increased sensitivity, specificity, and yield, while allowing assembly of reactions at room temperature. The use of Platinum® antibodies helps reduce PCR optimization requirements, reaction set-up and handling time, and contamination risk, thereby improving PCR results for templates up to 5 kb.

**Platinum® *Taq* DNA Polymerase High Fidelity** is a mixture of recombinant *Taq* DNA polymerase, *Pyrococcus* species *GB-D* polymerase, and Platinum® *Taq* antibody. Platinum® antibody complexes with *Taq* DNA polymerase and inhibits activity at ambient temperatures, allowing room-temperature setup. Activity is restored after the PCR denaturation step at 94°C, providing an automatic “hot start” for the enzyme and increasing specificity, sensitivity, and yield.

*Pyrococcus* species *GB-D* polymerase is a proofreading enzyme that possesses a 3’ to 5’ exonuclease activity. The enzyme mixture in Platinum® *Taq* DNA Polymerase High Fidelity results in a six-fold increase in fidelity over *Taq* DNA polymerase alone and allows amplification of simple and complex DNA templates over a large range of target sizes, up to 12 kb with no optimization.

# PCR and Quantitative PCR, continued



## Note

Since PCR is a powerful technique capable of amplifying trace amounts of DNA, take all appropriate precautions to avoid sample contamination.

Annealing and extension conditions are dependent on primer  $T_m$ , and should be determined independently for each reaction.

If PCR efficiency is not optimal, repeat the reaction with a primer titration from 100 to 500 nM (final conc.) in 100-nM increments.

## PCR — Targets Up to 1 kb

The following protocol uses Platinum® *Taq* DNA Polymerase in a standard PCR reaction. Adjust the reaction size as needed. Optimal reaction conditions—including incubation times and temperatures, and concentrations of enzyme, primers, and  $MgCl_2$ —may vary.

**Note:** A concentration of 1.5 mM  $MgCl_2$  is sufficient for most targets. For further optimization, prepare a titration from 1.5 mM to 3 mM in 0.25-mM increments.

1. Add the following components to a sterile 0.2- or 0.5-ml PCR tube or plate well at room temperature or on ice. For multiple reactions, prepare a master mix of common components.

Components	Volume	Final Conc.
10X PCR Buffer, Minus Mg	5 $\mu$ l	1X
10 mM dNTP mixture	1 $\mu$ l	0.2 mM each
50-mM $MgCl_2$	1.5 $\mu$ l	1.5 mM
Sense primer (10 $\mu$ M)	1 $\mu$ l	0.2 $\mu$ M
Antisense primer (10 $\mu$ M)	1 $\mu$ l	0.2 $\mu$ M
cDNA from Step 10, page 5	2 $\mu$ l	—
Platinum® <i>Taq</i> DNA Polymerase	0.4 $\mu$ l	2.0 units*
Autoclaved, distilled water	to 50 $\mu$ l	n/a

\*2.0 units are recommended for amplifying cDNA from the CellsDirect kit. In some cases, more enzyme may be required (up to 2.5 units).

2. Mix contents of the tubes and overlay with 50  $\mu$ l of mineral or silicone oil, if necessary.
3. Cap the tubes and centrifuge briefly to collect the contents.
4. Incubate tubes in a thermal cycler at 94°C for 30 seconds to 2 minutes to denature the template and activate the enzyme.
5. Perform 30–40 cycles of PCR amplification as follows:

Denature	94°C for 15–30 seconds
Anneal	55–65°C for 30 seconds
Extend	72°C for 1 minute per kb
6. Maintain the reaction at 4°C after cycling. The samples can be stored at -20°C until use.
7. Analyze the products by agarose gel electrophoresis and visualize by ethidium bromide staining. Use appropriate molecular weight standards.

Continued on next page

# PCR and Quantitative PCR, continued

## PCR — Targets Above 1 kb

The following protocol uses Platinum® *Taq* DNA Polymerase High Fidelity in a standard PCR reaction. Adjust the reaction size as needed. Optimal reaction conditions—including incubation times and temperatures, and the concentrations of Platinum® *Taq* DNA Polymerase High Fidelity, primers, MgSO<sub>4</sub>, and template DNA—may vary.

**Note:** A concentration of 2 mM MgSO<sub>4</sub> is sufficient for most targets. For further optimization, prepare a titration from 2 mM to 4 mM in 0.25-mM increments.

1. Add the following components to a sterile 0.2- or 0.5-ml PCR tube or plate well at room temperature or on ice. For multiple reactions, prepare a master mix of common components.

Component	Volume	Final Conc
10X High Fidelity PCR Buffer	5 µl	1X
10-mM dNTP mixture	1 µl	0.2 mM each
50-mM MgSO <sub>4</sub>	2 µl	2 mM
Sense primer (10 µM)	1 µl	0.2 µM
Antisense primer (10 µM)	1 µl	0.2 µM
cDNA from Step 10, page 5	≥1 µl	—
Platinum® <i>Taq</i> High Fidelity	0.2 µl	1.0 unit*
Autoclaved, distilled water	to 50 µl	Not applicable

\*1.0 unit is sufficient for amplifying most targets. In some cases, more enzyme may be required (up to 2.5 units).

2. Mix contents of the tubes and overlay with 50 µl of mineral or silicone oil, if necessary.
3. Cap the tubes and centrifuge briefly to collect the contents.
4. Incubate tubes in a thermal cycler at 94°C for 30 seconds to 2 minutes to denature the template and activate the enzyme.
5. Perform 30–40 cycles of PCR amplification as follows:

Denature

Anneal

Extend

94°C for 15–30 seconds

55–65°C for 30 seconds

68°C for 1 minute per kb
6. Maintain the reaction at 4°C after cycling. The samples can be stored at -20°C until use.
7. Analyze the products by agarose gel electrophoresis and visualize by ethidium bromide staining. Use appropriate molecular weight standards.

## Real-Time qPCR

We recommend Platinum® Quantitative PCR SuperMix-UDG for optimal results in real-time qPCR using LUX™ Primers or TaqMan® probes. We recommend Platinum® SYBR® Green qPCR SuperMix-UDG for optimal results using SYBR® Green I binding dye. See these product manuals for detailed protocols using first-strand cDNA in real-time qPCR on different instruments.

Continued on next page

# PCR and Quantitative PCR, continued

---

## PCR — Control Reaction

The following protocol uses Platinum® *Taq* DNA Polymerase and the Control Primers provided in the kit in a standard PCR reaction.

1. Prepare a PCR mixture for each control reaction from Step 10, page 5. For each control reaction, add the following to a sterile 0.2- or 0.5-ml PCR tube or plate well at either room temperature or on ice:

<u>Component</u>	<u>Volume</u>
DEPC-treated water	38.1 µl
10X PCR buffer minus Mg <sup>++</sup>	5 µl
50 mM MgCl <sub>2</sub>	1.5 µl
10 mM dNTP mix	1 µl
Forward Control Primer (10 µM)	1 µl
Reverse Control Primer (10 µM)	1 µl
cDNA from control RNA/ negative RT control, Step 10, page 5	2 µl
Platinum® <i>Taq</i> DNA polymerase (5 units/µl)	<u>0.4 µl</u>
final volume	50 µl

2. Mix the contents of the tube. Centrifuge briefly to collect the reaction components.
3. Place reaction mixture in preheated (94°C) thermal cycler. Perform an initial denaturation step: 94°C for 2 min.
4. Perform 40 cycles of PCR:

Denature

94°C for 15 sec

Anneal

60°C for 30 sec

Extend

72°C for 1 min

**Note:** For slow-ramping thermal cyclers, follow manufacturer’s directions.

5. Upon completion, maintain reactions at 4°C.
  6. Analyze 10 µl of each sample using agarose gel electrophoresis and ethidium bromide staining. A 1.18-kb band corresponding to at least 25 ng of product should be visible for the control sample. No band should be visible for the negative RT control sample.
-

# Troubleshooting

Problem	Possible Cause	Suggested Solution
No bands after electrophoretic analysis of amplified products	Procedural error	Confirm that all steps were followed. Use the Control RNA to verify the efficiency of the first-strand reaction (see the next page on troubleshooting with the Control RNA).
	RNA is degraded	<p>Add control total HeLa RNA to sample to determine if RNase is present in the first-strand reaction.</p> <p>Confirm that RNaseOUT™ was added at the appropriate steps in the protocol.</p> <p>A longer DNase I digestion can hydrolyze the RNA in the sample. Use a digestion time of &lt;10 minutes.</p> <p>Maintain aseptic conditions to prevent RNase contamination.</p>
	Target mRNA contains strong transcriptional pauses	<p>Use random hexamers (Cat. no. 48190-011) instead of oligo(dT)<sub>20</sub> in the first-strand reaction.</p> <p>Maintain an elevated temperature after the annealing step.</p> <p>Increase the temperature of first-strand reaction (up to 55°C).</p> <p>Use PCR primers closer to the 3' terminus of the target cDNA.</p>
	Too much first-strand product was used in PCR	Use no more than 5 µl of the first-strand product in PCR.
Unexpected bands after electrophoretic analysis	Contamination by genomic DNA	<p>Do not omit the DNase Digestion step on page 4. For larger samples (&gt;1,000 cells), use a longer DNase I incubation time, i.e., up to 10 minutes.</p> <p>Design primers that anneal to sequence in exons on both sides of an intron or exon/exon boundary of the mRNA to of amplified allow differentiation between amplification of cDNA and products potential contaminating genomic DNA.</p> <p>To test if products were derived from DNA, prepare a negative RT control.</p>
	Nonspecific annealing of primers	<p>Vary the annealing conditions. Use Platinum® Taq DNA Polymerase for automatic hot-start PCR.</p> <p>Optimize magnesium concentration for each template and primer combination.</p>



## References

---

- Berger, S.L. and Kimmel, A.R. (1987) *Methods Enzymol* 152, 316.
- Bracete, A.M., Mertz, L.M., Fox, D.K. (1999) *Focus*® 21, 38.
- Chomczynski, P. (1993) *Biotechniques* Vol. 15, 532.
- Chomczynski, P. and Sacchi, N. (1987) *Anal. Biochem.* 162, 156.
- Compton, T. (1990) in *PCR Protocols: A Guide to Methods and Applications* (Innis, M., Gelfand, D., Sninsky, J., and White, T., eds.), p. 39, Academic Press, Inc.
- D'Alessio, J. M., Gruber, C. E., Cain, C., and Noon, M. C. (1990) *Focus*® 12, 47.
- Frohman, M.A., Dush, M.K, and Martin, G.R. (1988) *Proc. Nat. Acad. Sci USA* 85, 8998.
- Gerard, G.F. (1994) *Focus*® 16, 102.
- Gerard, G.F., D'Alessio, J.M., and Kotewicz, M.L. (1989) *Focus*® 11, 66.
- Gerard, G.F., Schmidt, B.J., Kotewicz, M.L., and Campbell, J.H. (1992) *Focus*® 14, 91.
- Hu, A.W., D'Alessio, J.M., Gerard, G.F., and Kullman, J. (1991) *Focus*® 13, 26.
- Lee, C.C. and Caskey, T. (1990) in *PCR Protocols: A Guide to Methods and Applications* (Innis, M., Gelfand, D., Sninsky, J., and White, T., eds.), p. 46, Academic Press, Inc.
- Sambrook J., Fritsch, E.F., and Maniatis, T. (1989) *Molecular Cloning: A Laboratory Manual*, Cold Spring Harbor Laboratory Press, Cold Spring Harbor.
- Simms, D., Cizdziel, P.E., and Chomczynski, P. (1993) *Focus*® 15, 99.
- Westfall, B., Sitaraman, K., Solus, J., Hughes, J., and Rashtchian, A. (1997) *Focus*® 19, 46.
- Westfall, B., Sitaraman, K., Berninger, M., and Mertz, L.M. (1995) *Focus*® 17, 62.
- Westfall, B., Sitaraman, K., Lee, J., Borman, J. and Rashtchian, A. (1999) *Focus*® 21, 49.
- Takagi, M., Nishioka, M., Kakiyama, H., Kitabayashi, M., Inoue, H., Kawakami, B., Oka, M., and Imanaka, T. (1997) *Appl. Environ. Microbiol.* 63, 4504.
- Sitaraman, K., Darfler, M., and Westfall, B. (1999) *Focus*® 21, 10.
- Nathan, M., Mertz, L., Fox, D. (1995) *Focus*® 17, 78.
- Schwabe, W., Lee, J.E., Nathan, M., Xu, R.H., Sitaraman, K., Smith, M., Potter, R.J., Rosenthal, K., Rashtchian, A., Gerard, G.F. (1998) *Focus*® 20, 30.
-

## Related Products

Product	Size	Cat. No.
<b><u>Products for Analysis:</u></b>		
E-Gel® Pre-cast Agarose Gels		
0.8% Starter Pak	9 gels and base	G5000-08
1.2% Starter Pak	9 gels and base	G5000-01
2% Starter Pak	9 gels and base	G5000-02
4% Starter Pak	9 gels and base	G5000-04
UltraPure™ Agarose	100 g	15510-019
	500 g	15510-027
UltraPure™ Agarose 1000	100 g	10975-035
100-bp DNA Ladder	50 µg	15628-019
123-bp DNA Ladder	100 µg	15613-011
	250 µg	15613-029
1-Kb Plus DNA Ladder	250 µg	10787-018
	1,000 µg	10787-026
<b><u>Enzymes for Amplification:</u></b>		
Platinum® <i>Taq</i> DNA Polymerase	100 reactions	10966-018
	250 reactions	10966-026
	500 reactions	10966-034
	5000 reactions	10966-083
Platinum® <i>Taq</i> DNA Polymerase High Fidelity	100 reactions	11304-011
	500 reactions	11304-029
	5,000 reactions	11304-102
Platinum® Quantitative PCR	100 rxns	11730-017
SuperMix-UDG	250 reactions	11730-025
PCR SuperMix High Fidelity	100 reactions	10790-020
Platinum® PCR SuperMix	100 reactions	11306-016

## Purchaser Notification

---

**Limited Use  
Label License  
No. 4: Products  
for PCR that  
include no  
rights to  
perform PCR**

This product is optimized for use in the Polymerase Chain Reaction (PCR) covered by patents owned by Roche Molecular Systems, Inc. and F. Hoffmann-La Roche, Ltd. ("Roche"). No license under these patents to use the PCR process is conveyed expressly or by implication to the purchaser by the purchase of this product. A license to use the PCR process for certain research and development activities accompanies the purchase of certain reagents from licensed suppliers such as Invitrogen, when used in conjunction with an Authorized Thermal Cycler, or is available from Applied Biosystems. Further information on purchasing licenses to practice the PCR process may be obtained by contacting the Director of Licensing at Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404 or at Roche Molecular Systems, Inc., 1145 Atlantic Avenue, Alameda, California 94501.

---

**Limited Use  
Label License  
No. 138:  
SuperScript™ III  
Reverse  
Transcriptase**

The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) not to transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. Invitrogen Corporation will not assert a claim against the buyer of infringement of patents owned by Invitrogen Corporation and claiming this product based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of this limited use statement, Invitrogen is willing to accept return of the product with a full refund. For information on purchasing a license to this product for purposes other than research, contact Licensing Department, Invitrogen Corporation, 1600 Faraday Avenue, Carlsbad, California 92008. Phone (760) 603-7200. Fax (760) 602-6500.

---

*Continued on next page*

## Purchaser Notification, Continued

---

### **Limited Use Label License No: 18: RNaseOUT™ Ribonuclease Inhibitor**

This product is the subject of U.S. Patent No. 5,965,399 owned by Invitrogen Corporation. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) to not transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. Invitrogen Corporation will not assert a claim against the buyer of infringement of the above patents based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of this limited use statement, Invitrogen is willing to accept return of the product with a full refund. For information on purchasing a license to this product for purposes other than research, contact Licensing Department, Invitrogen Corporation, 1600 Faraday Avenue, Carlsbad, California 92008. Phone (760) 603-7200. Fax (760) 602-6500.

---

### **Trademarks**

TaqMan® is a registered trademark of Roche Molecular Systems, Inc.

---

## Technical Service

---

### World Wide Web



Visit the Invitrogen Web Resource using your World Wide Web browser. At the site, you can:

- Get the scoop on our hot new products and special product offers
- View and download vector maps and sequences
- Download manuals in Adobe® Acrobat® (PDF) format
- Explore our catalog with full color graphics
- Obtain citations for Invitrogen products
- Request catalog and product literature

Once connected to the Internet, launch your Web browser (Internet Explorer 5.0 or newer or Netscape 4.0 or newer), then enter the following location (or URL):

<http://www.invitrogen.com>

...and the program will connect directly. Click on underlined text or outlined graphics to explore. Don't forget to put a bookmark at our site for easy reference!

---

### Contact Us

For more information or technical assistance, call, write, fax, or email. Additional international offices are listed on our Web page ([www.invitrogen.com](http://www.invitrogen.com)).

#### Corporate Headquarters:

Invitrogen Corporation  
1600 Faraday Avenue  
Carlsbad, CA 92008 USA  
Tel: 1 760 603 7200  
Tel (Toll Free): 1 800 955 6288  
Fax: 1 760 602 6500  
E-mail:  
[tech\\_service@invitrogen.com](mailto:tech_service@invitrogen.com)

#### Japanese Headquarters:

Invitrogen Japan K.K.  
Nihonbashi Hama-Cho Park  
Bldg. 4F  
2-35-4, Hama-Cho,  
Nihonbashi  
Tel: 81 3 3663 7972  
Fax: 81 3 3663 8242  
E-mail:  
[jpinfo@invitrogen.com](mailto:jpinfo@invitrogen.com)

#### European Headquarters:

Invitrogen Ltd  
Inchinnan Business Park  
3 Fountain Drive  
Paisley PA4 9RF, UK  
Tel: +44 (0) 141 814 6100  
Tech Fax: +44 (0) 141 814 6117  
E-mail:  
[eurotech@invitrogen.com](mailto:eurotech@invitrogen.com)

---

### MSDS Requests

To request an MSDS, visit our Web site at [www.invitrogen.com](http://www.invitrogen.com). On the home page, go to 'Technical Resources', select 'MSDS', and follow instructions on the page.

---

*Continued on next page*

## Technical Service, Continued

---

### Limited Warranty

Invitrogen is committed to providing our customers with high-quality goods and services. Our goal is to ensure that every customer is 100% satisfied with our products and our service. If you should have any questions or concerns about an Invitrogen product or service, contact our Technical Service Representatives.

Invitrogen warrants that all of its products will perform according to specifications stated on the certificate of analysis. The company will replace, free of charge, any product that does not meet those specifications. This warranty limits Invitrogen Corporation's liability only to the cost of the product. No warranty is granted for products beyond their listed expiration date. No warranty is applicable unless all product components are stored in accordance with instructions. Invitrogen reserves the right to select the method(s) used to analyze a product unless Invitrogen agrees to a specified method in writing prior to acceptance of the order.

Invitrogen makes every effort to ensure the accuracy of its publications, but realizes that the occasional typographical or other error is inevitable. Therefore Invitrogen makes no warranty of any kind regarding the contents of any publications or documentation. If you discover an error in any of our publications, please report it to our Technical Service Representatives.

**Invitrogen assumes no responsibility or liability for any special, incidental, indirect or consequential loss or damage whatsoever. The above limited warranty is sole and exclusive. No other warranty is made, whether expressed or implied, including any warranty of merchantability or fitness for a particular purpose.**

---

©2004 Invitrogen Corporation. All rights reserved.

For research use only. Not intended for any animal or human therapeutic or diagnostic use.





**Corporate Headquarters:**

Invitrogen Corporation  
1600 Faraday Avenue  
Carlsbad, California 92008  
Tel: 1 760 603 7200  
Tel (Toll Free): 1 800 955 6288  
Fax: 1 760 602 6500  
Email: [tech\\_service@invitrogen.com](mailto:tech_service@invitrogen.com)

**European Headquarters:**

Invitrogen Ltd  
Inchinnan Business Park  
3 Fountain Drive  
Paisley PA4 9RF, UK  
Tel (Free Phone Orders): 0800 269 210  
Tel (General Enquiries): +44 (0) 141 814 6100  
Fax: +44 (0) 141 814 6260  
Email: [eurotech@invitrogen.com](mailto:eurotech@invitrogen.com)

**International Offices:**

Argentina 5411 4556 0844  
Australia 1 800 331 627  
Brazil 55 11 5051 7422  
Canada 800 263 6236  
China 10 6849 2578  
Hong Kong 2407 8450  
Japan 03 3663 7974  
New Zealand 0800 600 200  
Singapore 65 686 186 38  
Taiwan 2 2651 6156

For other countries see our website

[www.invitrogen.com](http://www.invitrogen.com)

